

# **KMS Deep Vein Thrombosis (DVT) Pump**



- Intended for prophylaxis of deep vein thrombosis
- Enhances faster healing & tissue regeneration
- Weighs less than 5 lbs.
- Garments one size fits all & disposable
- For use in facilities and at home



**CUSTOMER SERVICE**

**1-800-488-6535 ext 2 • Fax 310-533-8810**

**[www.kingmedicalsupply.com](http://www.kingmedicalsupply.com)**

# KMS DVT PUMP

Intended Use	DVT Prophylaxis	
Model	IC - 1545- DL	
Garments	DVT Knee-High Bilateral Garment single patient use	
Timing	ON	OFF
	13-17 seconds	43-47 seconds
Pressure	50mmHg preset override available	
Alarm Pressure	10-20mmHg	

## CONTRAINDICATIONS

Use of this device is contraindicated for patients with the following condition:

- Any local condition of the extremity that would interfere with its application, including but not limited to:
  - Dermatitis
  - Immediately following vein ligation
  - Gangrene
  - Skin grafts
  - Casts or splints
- Severe arteriosclerosis or other ischemic vascular disease
- Pulmonary edema
- Congestive heart failure
- Suspicion or confirmation of the presence of existing DVT
- Infections of the limb, including cellulitis without appropriate antibiotic coverage for at least 72 hours prior to treatment
- Presence of lymphangiosarcoma
- Inflammatory phlebitis or episodes of pulmonary embolism



King Medical Supply  
[www.kingmedicalsupply.com](http://www.kingmedicalsupply.com)  
 800-488-6535