





KMS FULL ELECTRIC BED SYSTEM



STANDARD FEATURES

- Fully electric adjustments for bed height, upper body positions, and knee positions all controlled with a hand pendant.
- Welded angle steel construction with reinforced corner plates provide maximum stability.
- U-Roll steel bed surface provides maximum support with no mattress sagging.
- Two piece split frame design is easy to use, set up and operate.
- Emergency hand crank to manually raise and lower any parts of bed.
- 3" swivel, locking casters provide increased security.
- Accommodates weights up to 460lbs.
- Bariatric Sizes Available

| FABRIC HEIGHT | OVERALL WIDTH | OVERALL LENGTH | MATTRESS SIZE |
|-------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|
|  |  |  |  |
| 16" to 24" | 36" | 88" | 35" X 80" |

1-800-488-6535



www.kingmedicalsupply.com